

Recombinant Human SENP8

Catalog No: #GP12001



Package Size: #GP12001-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SENP8
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	SUMO/sentrin specific peptidase family member 8
Other Names	DEN1; NEDP1; PRSC2
Accession No.	Swissprot:Q96LD8Gene Accession:BC031411
Uniprot	Q96LD8
GeneID	123228;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a cysteine protease that is a member of the sentrin-specific protease family. The encoded protein is involved in processing and deconjugation of the ubiquitin-like protein termed, neural precursor cell expressed developmentally downregulated 8. Alternate splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only