

Recombinant Human SENP2

Catalog No: #GP12003



Package Size: #GP12003-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SENP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 189 amino acids of human SUMO1/sentrin/SMT3 specific peptidase 2
Target Name	SUMO1/sentrin/SMT3 specific peptidase 2
Other Names	AXAM2; SMT3IP2
Accession No.	Swissprot:Q9HC62Gene Accession:BC040609
Uniprot	Q9HC62
GeneID	59343;
Storage	-20~-80°C, pH 7.6 PBS

Background

SUMO1 (UBL1; MIM 601912) is a small ubiquitin-like protein that can be covalently conjugated to other proteins. SENP2 is one of a group of enzymes that process newly synthesized SUMO1 into the conjugatable form and catalyze the deconjugation of SUMO1-containing species.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only