

Recombinant Human SENP7

Catalog No: #GP12004



Package Size: #GP12004-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SENP7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 801-1000 amino acids of human SUMO1/sentrin specific peptidase 7
Target Name	SUMO1/sentrin specific peptidase 7
Accession No.	Swissprot:Q9BQF6Gene Accession:BC129988
Uniprot	Q9BQF6
GeneID	57337;
Storage	-20~-80°C, pH 7.6 PBS

Background

The reversible posttranslational modification of proteins by the addition of small ubiquitin-like SUMO proteins (see SUMO1; MIM 601912) is required for many cellular processes. SUMO-specific proteases, such as SENP7, process SUMO precursors to generate a C-terminal diglycine motif required for the conjugation reaction.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only