

Recombinant Human SERPINB4

Catalog No: #GP12006



Package Size: #GP12006-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SERPINB4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-390 amino acids of human SERPINB4
Target Name	serpin family B member 4
Other Names	PI11; SCCA1; SCCA2; LEUPIN; SCCA-2
Accession No.	Swissprot:P48594Gene Accession:BC017401
Uniprot	P48594
GeneID	6318;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a member of the serpin family of serine protease inhibitors. The encoded protein is highly expressed in many tumor cells and can inactivate granzyme M, an enzyme that kills tumor cells. This protein, along with serpin B3, can be processed into smaller fragments that aggregate to form an autoantigen in psoriasis, probably by causing chronic inflammation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only