

# Recombinant Human SH3KBP1

Catalog No: #GP12014



Package Size: #GP12014-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SH3KBP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 215 amino acids of human SH3-domain kinase binding protein 1
Target Name	SH3-domain kinase binding protein 1
Other Names	HSB1; CIN85; GIG10; HSB-1; MIG18; CD2BP3
Accession No.	Swissprot:Q96B97Gene Accession:BC015806
Uniprot	Q96B97
GeneID	30011;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes an adapter protein that contains three N-terminal Src homology domains, a proline rich region and a C-terminal coiled-coil domain. The encoded protein facilitates protein-protein interactions and has been implicated in numerous cellular processes including apoptosis, cytoskeletal rearrangement, cell adhesion and in the regulation of clathrin-dependent endocytosis. Alternate splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only