

# Recombinant Human SIAH1

Catalog No: #GP12017



Package Size: #GP12017-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SIAH1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	siah E3 ubiquitin protein ligase 1
Other Names	SIAH1A
Accession No.	Swissprot:Q8IUQ4Gene Accession:BC035562
Uniprot	Q8IUQ4
GeneID	6477;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a protein that is a member of the seven in absentia homolog (SIAH) family. The protein is an E3 ligase and is involved in ubiquitination and proteasome-mediated degradation of specific proteins. The activity of this ubiquitin ligase has been implicated in the development of certain forms of Parkinson's disease, the regulation of the cellular response to hypoxia and induction of apoptosis. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only