

Recombinant Human SRPRB

Catalog No: #GP12021



Package Size: #GP12021-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SRPRB
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	signal recognition particle receptor, B subunit
Other Names	APMCF1
Accession No.	Swissprot:Q9Y5M8Gene Accession:BC065299
Uniprot	Q9Y5M8
GeneID	58477;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene has similarity to mouse protein which is a subunit of the signal recognition particle receptor (SR). This subunit is a transmembrane GTPase belonging to the GTPase superfamily. It anchors alpha subunit, a peripheral membrane GTPase, to the ER membrane. SR is required for the cotranslational targeting of both secretory and membrane proteins to the ER membrane.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only