

Recombinant Human SLC25A26

Catalog No: #GP12026



Package Size: #GP12026-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SLC25A26
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 186 amino acids of human solute carrier family 25 (S-adenosylmethionine carrier), member 26
Target Name	solute carrier family 25 (S-adenosylmethionine carrier), member 26
Other Names	SAMC
Accession No.	Swissprot:Q70HW3Gene Accession:BC012852
Uniprot	Q70HW3
GenID	115286;
Storage	-20~-80°C, pH 7.6 PBS

Background

Mitochondrial carriers, including SLC25A26, are a family of transport proteins found mostly in the inner membranes of mitochondria. They shuttle metabolites and cofactors through the mitochondrial membrane?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only