

Recombinant Human SLITRK4

Catalog No: #GP12032



Package Size: #GP12032-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SLITRK4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 197 amino acids of human SLIT and NTRK-like family, member 4
Target Name	SLIT and NTRK-like family, member 4
Accession No.	Swissprot:Q8IW52Gene Accession:BC040986
Uniprot	Q8IW52
GeneID	139065;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a transmembrane protein belonging to the the SLITRK family. These family members include two N-terminal leucine-rich repeat domains similar to those found in the axonal growth-controlling protein SLIT, as well as C-terminal regions similar to neurotrophin receptors. Studies of an homologous protein in mouse suggest that this family member functions to suppress neurite outgrowth. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only