

Recombinant Human TAGLN

Catalog No: #GP12035



Package Size: #GP12035-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TAGLN
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	transgelin
Other Names	SM22; SMCC; TAGLN1; WS3-10
Accession No.	Swissprot:Q01995Gene Accession:BC004927
Uniprot	Q01995
GeneID	6876;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is a transformation and shape-change sensitive actin cross-linking/gelling protein found in fibroblasts and smooth muscle. Its expression is down-regulated in many cell lines, and this down-regulation may be an early and sensitive marker for the onset of transformation. A functional role of this protein is unclear. Two transcript variants encoding the same protein have been found for this gene.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only