

# Recombinant Human HELLS

Catalog No: #GP12036



Package Size: #GP12036-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human HELLS
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 591-838 amino acids of human HELLS
Target Name	helicase, lymphoid-specific
Other Names	LSH; PASG; SMARCA6; Nbla10143
Accession No.	Swissprot:Q9NRZ9Gene Accession:BC030963
Uniprot	Q9NRZ9
GeneID	3070;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only