

# Recombinant Human SMYD3

Catalog No: #GP12038



Package Size: #GP12038-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human SMYD3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 228 amino acids of human SET and MYND domain containing 3
Target Name	SET and MYND domain containing 3
Other Names	KMT3E; ZMYND1; ZNFN3A1; bA74P14.1
Accession No.	Swissprot:Q9H7B4Gene Accession:BC031010
Uniprot	Q9H7B4
GeneID	64754;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a histone methyltransferase which functions in RNA polymerase II complexes by an interaction with a specific RNA helicase. Multiple transcript variants encoding different isoforms have been found for this gene.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only