

Recombinant Human SMYD4

Catalog No: #GP12039



Package Size: #GP12039-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SMYD4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 204 amino acids of human SMYD4
Target Name	SET and MYND domain containing 4
Other Names	ZMYND21
Accession No.	Swissprot:Q8IYR2Gene Accession:BC035077
Uniprot	Q8IYR2
GeneID	114826;
Storage	-20~-80°C, pH 7.6 PBS

Background

The specific function of the protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only