

Recombinant Human SNAI3

Catalog No: #GP12040



Package Size: #GP12040-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SNAI3
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	snail family zinc finger 3
Other Names	SMUC; SNAIL3; ZNF293; Zfp293
Accession No.	Swissprot:Q3KNW1 Gene Accession:BC107058
Uniprot	Q3KNW1
GeneID	333929;
Storage	-20~-80°C, pH 7.6 PBS

Background

SNAI3 is a member of the SNAIL gene family, named for the *Drosophila* snail gene, which plays roles in mesodermal formation during embryogenesis

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only