

Recombinant Human SNAP23

Catalog No: #GP12041



Package Size: #GP12041-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SNAP23
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	synaptosomal-associated protein, 23kDa
Other Names	SNAP-23; SNAP23A; SNAP23B; HsT17016
Accession No.	Swissprot:O00161Gene Accession:BC003686
Uniprot	O00161
GeneID	8773;
Storage	-20~-80°C, pH 7.6 PBS

Background

Specificity of vesicular transport is regulated, in part, by the interaction of a vesicle-associated membrane protein termed synaptobrevin/VAMP with a target compartment membrane protein termed syntaxin. These proteins, together with SNAP25 (synaptosome-associated protein of 25 kDa), form a complex which serves as a binding site for the general membrane fusion machinery.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only