

Recombinant Human SNX20

Catalog No: #GP12052



Package Size: #GP12052-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human SNX20 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to N terminal 95 amino acids of human sorting nexin 20 |
| Target Name | sorting nexin 20 |
| Other Names | SLIC1 |
| Accession No. | Swissprot:Q7Z614Gene Accession:BC063423 |
| Uniprot | Q7Z614 |
| GeneID | 124460; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

SNX20 interacts with the cytoplasmic domain of PSGL1 (SELPLG; MIM 600738) and cycles PSGL1 into endosomes.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only