

Recombinant Human SPRY4

Catalog No: #GP12054



Package Size: #GP12054-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SPRY4
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	sprouty RTK signaling antagonist 4
Other Names	HH17
Accession No.	Swissprot:Q9C004Gene Accession:BC125095
Uniprot	Q9C004
GeneID	81848;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of a family of cysteine- and proline-rich proteins. The encoded protein is an inhibitor of the receptor-transduced mitogen-activated protein kinase (MAPK) signaling pathway. Activity of this protein impairs the formation of active GTP-RAS. Nucleotide variation in this gene has been associated with hypogonadotropic hypogonadism 17 with or without anosmia. Alternative splicing results in a multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only