

Recombinant Human SSX1

Catalog No: #GP12060



Package Size: #GP12060-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SSX1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	synovial sarcoma, X breakpoint 1
Other Names	SSRC; CT5.1
Accession No.	Swissprot:Q16384?Gene Accession:BC128611
Uniprot	Q16384
GeneID	6756;
Storage	-20~-80°C, pH 7.6 PBS

Background

The product of this gene belongs to the family of highly homologous synovial sarcoma X (SSX) breakpoint proteins. These proteins may function as transcriptional repressors. They are also capable of eliciting spontaneous humoral and cellular immune responses in cancer patients, and are potentially useful targets in cancer vaccine-based immunotherapy.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only