

Recombinant Human SDCBP

Catalog No: #GP12062



Package Size: #GP12062-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SDCBP
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	syndecan binding protein (syntenin)
Other Names	ST1; SYCL; MDA-9; TACIP18
Accession No.	Swissprot:O00560Gene Accession:BC113674
Uniprot	O00560
GeneID	6386;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene was initially identified as a molecule linking syndecan-mediated signaling to the cytoskeleton. The syntenin protein contains tandemly repeated PDZ domains that bind the cytoplasmic, C-terminal domains of a variety of transmembrane proteins. This protein may also affect cytoskeletal-membrane organization, cell adhesion, protein trafficking, and the activation of transcription factors.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only