

Recombinant Human MED22

Catalog No: #GP12075



Package Size: #GP12075-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human MED22 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Full length fusion protein |
| Target Name | mediator complex subunit 22 |
| Other Names | MED24; SURF5; surf-5 |
| Accession No. | Swissprot:Q15528Gene Accession:BC040111 |
| Uniprot | Q15528 |
| GeneID | 6837; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes a protein component of the mediator complex, which functions in the regulation of transcription by bridging interactions between gene-specific regulatory factors, RNA polymerase II, and general transcription factors. Alternatively spliced transcript variants encoding different isoforms have been observed.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only