

## Recombinant Human ARMC10

Catalog No: #GP12076



Package Size: #GP12076-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ARMC10
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	armadillo repeat containing 10
Other Names	SVH; PNAS112; PNAS-112; PSEC0198
Accession No.	Swissprot:Q8N2F6Gene Accession:BC003586
Uniprot	Q8N2F6
GeneID	83787;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a protein that contains an armadillo repeat and transmembrane domain. The encoded protein decreases the transcriptional activity of the tumor suppressor protein p53 through direct interaction with the DNA-binding domain of p53, and may play a role in cell growth and survival. Upregulation of this gene may play a role in hepatocellular carcinoma. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 3.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only