

Recombinant Human SYTL5

Catalog No: #GP12077



Package Size: #GP12077-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SYTL5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 230 amino acids of human synaptotagmin-like 5
Target Name	synaptotagmin-like 5
Other Names	slp5
Accession No.	Swissprot:Q8TDW5Gene Accession:BC131585
Uniprot	Q8TDW5
GeneID	94122;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene belongs to the synaptotagmin-like (Slp) protein family, which contains a unique homology domain at the N-terminus, referred to as the Slp homology domain (SHD). The SHD functions as a binding site for Rab27A, which plays a role in protein transport. Expression of this gene is restricted to placenta and liver, suggesting that it might be involved in Rab27A-dependent membrane trafficking in specific tissues.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only