

# Recombinant Human STX7

Catalog No: #GP12078



Package Size: #GP12078-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human STX7
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 131 amino acids of human STX7
Target Name	syntaxin 7
Accession No.	Swissprot:O15400Gene Accession:BC011975
Uniprot	O15400
GeneID	8417;
Storage	-20~-80°C, pH 7.6 PBS

## Background

May be involved in protein trafficking from the plasma membrane to the early endosome (EE) as well as in homotypic fusion of endocytic organelles. Mediates the endocytic trafficking from early endosomes to late endosomes and lysosomes.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only