

## Recombinant Human STX8

Catalog No: #GP12079



Package Size: #GP12079-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human STX8
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 146 amino acids of human syntaxin 8
Target Name	syntaxin 8
Other Names	CARB
Accession No.	Swissprot:Q9UNK0Gene Accession:BC009713
Uniprot	Q9UNK0
GeneID	9482;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The gene is a member of the syntaxin family. The encoded protein is involved in protein trafficking from early to late endosomes via vesicle fusion and exocytosis. A related pseudogene has been identified on chromosome 12. Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only