

## Recombinant Human STX12

Catalog No: #GP12080



Package Size: #GP12080-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human STX12
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	syntaxin 12
Other Names	STX13; STX14
Accession No.	Swissprot:Q86Y82Gene Accession:BC046999
Uniprot	Q86Y82
GeneID	23673;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Syntaxin 12, also known as STX12, STX13 or STX14, is a 276 amino acid single-pass membrane protein that contains one t-SNARE coiled-coil homology domain and belongs to the Syntaxin family. Syntaxin 12 regulates protein transport between late endosomes and the trans-Golgi network and interacts with ABC1 (ATP-binding cassette transporter A1), a protein that facilitates cellular release of choline-phospholipids and cholesterol to apoA-I (apolipoprotein A-I).?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only