

Recombinant Human STX16

Catalog No: #GP12081



Package Size: #GP12081-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human STX16
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	syntaxin 16
Other Names	SYN16
Accession No.	Swissprot:O14662Gene Accession:BC073876
Uniprot	O14662
GeneID	8675;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that is a member of the syntaxin or t-SNARE (target-SNAP receptor) family. These proteins are found on cell membranes and serve as the targets for V-SNARES (vesicle-SNAP receptors) permitting specific synaptic vesicle docking and fusion. A microdeletion in the region of chromosome 20 where this gene is located has been associated with pseudohypoparathyroidism type 1b.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only