

Recombinant Human STX11

Catalog No: #GP12085



Package Size: #GP12085-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human STX11
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	syntaxin 11
Other Names	FHL4; HLH4; HPLH4
Accession No.	Swissprot:O75558Gene Accession:BC033519
Uniprot	O75558
GeneID	8676;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the syntaxin family. Syntaxins have been implicated in the targeting and fusion of intracellular transport vesicles. This family member may regulate protein transport among late endosomes and the trans-Golgi network. Mutations in this gene have been associated with familial hemophagocytic lymphohistiocytosis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only