

Recombinant Human STX19

Catalog No: #GP12086



Package Size: #GP12086-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human STX19
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	syntaxin 19
Accession No.	Swissprot:Q8N4C7?Gene Accession:BC034696
Uniprot	Q8N4C7
GeneID	415117;
Storage	-20~-80°C, pH 7.6 PBS

Background

Syntaxin 19, also known as STX19, is a 294 amino acid peripheral membrane protein that contains one t-SNARE coiled-coil homology domain and belongs to the syntaxin family, suggesting a role in synaptic vesicle fusion. The gene encoding Syntaxin 19 maps to human chromosome 3, which houses over 1,100 genes, including a chemokine receptor (CKR) gene cluster and a variety of human cancer-related gene loci.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only