

Recombinant Human TBL1XR1

Catalog No: #GP12094



Package Size: #GP12094-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TBL1XR1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 214 amino acids of human transducin (beta)-like 1 X-linked receptor 1
Target Name	transducin (beta)-like 1 X-linked receptor 1
Other Names	C21; DC42; IRA1; TBLR1
Accession No.	Swissprot:Q9BZK7Gene Accession:BC113421
Uniprot	Q9BZK7
GenID	79718;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene has sequence similarity with members of the WD40 repeat-containing protein family. The WD40 group is a large family of proteins, which appear to have a regulatory function. It is believed that the WD40 repeats mediate protein-protein interactions and members of the family are involved in signal transduction, RNA processing, gene regulation, vesicular trafficking, cytoskeletal assembly and may play a role in the control of cytotypic differentiation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only