

Recombinant Human TBX18

Catalog No: #GP12095



Package Size: #GP12095-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TBX18
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 257 amino acids of human T-box 18
Target Name	T-box 18
Accession No.	Swissprot:O95935Gene Accession:BC132715
Uniprot	O95935
GeneID	9096;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene codes for a member of an evolutionarily conserved family of transcription factors that plays a crucial role in embryonic development. The family is characterized by the presence of the DNA-binding T-box domain and is divided into five sub-families based on sequence conservation in this domain. The encoded protein belongs to the vertebrate specific Tbx1 sub-family. The protein acts as a transcriptional repressor by antagonizing transcriptional activators in the T-box family. The protein forms homo- or heterodimers with other transcription factors of the T-box family or other transcription factors.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only