

# Recombinant Human TBX10

Catalog No: #GP12096



Package Size: #GP12096-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TBX10
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 185 amino acids of human T-box 10
Target Name	T-box 10
Other Names	TBX7; TBX13
Accession No.	Swissprot:O75333Gene Accession:BC113485
Uniprot	O75333
GeneID	347853;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the T-box family of transcription factors. These transcription factors share a DNA-binding domain called the T-box, and play a role in several developmental processes including early embryonic cell fate and organogenesis. The encoded protein is a member of the T-box 1 subfamily. Mutations in this gene are thought to be a cause of isolated cleft lip with or without cleft palate.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only