

Recombinant Human RASD2

Catalog No: #GP12099



Package Size: #GP12099-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RASD2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	RASD family, member 2
Other Names	Rhes; TEM2; MGC:4834
Accession No.	Swissprot:Q96D21 Gene Accession:BC013419
Uniprot	Q96D21
GeneID	23551;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a Ras-related protein that enriched in striatum. The product of this gene binds to GTP and possesses intrinsic GTPase activity. The gene belongs to the Ras superfamily of small GTPases. The exact function of this gene is unknown, but most striatum-specific mRNAs characterized to date encode components of signal transduction cascades.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only