

Recombinant Human TFDP2

Catalog No: #GP12100



Package Size: #GP12100-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TFDP2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	transcription factor Dp-2 (E2F dimerization partner 2)
Other Names	DP2
Accession No.	Swissprot:Q14188Gene Accession:BC021113
Uniprot	Q14188
GeneID	7029;
Storage	-20~-80°C, pH 7.6 PBS

Background

The gene is a member of the transcription factor DP family. The encoded protein forms heterodimers with the E2F transcription factors resulting in transcriptional activation of cell cycle regulated genes. Alternative splicing results in multiple transcript variants.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only