

# Recombinant Human THAP10

Catalog No: #GP12102



Package Size: #GP12102-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human THAP10
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	THAP domain containing 10
Accession No.	Swissprot:Q9P2Z0Gene Accession:BC027857
Uniprot	Q9P2Z0
GeneID	56906;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of a family of proteins sharing an N-terminal Thanatos-associated domain. The Thanatos-associated domain contains a zinc finger signature similar to DNA-binding domains. This gene is part of a bidirectional gene pair on the long arm of chromosome 15 that is regulated by estrogen and may play a role in breast cancer.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only