

# Recombinant Human THAP2

Catalog No: #GP12103



Package Size: #GP12103-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human THAP2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	THAP domain containing, apoptosis associated protein 2
Accession No.	Swissprot:Q9H0W7Gene Accession:BC008358
Uniprot	Q9H0W7
GeneID	83591;
Storage	-20~-80°C, pH 7.6 PBS

## Background

THAP2, also named as THAP domain-containing protein 2, is a 228 amino acid protein, which contains one THAP-type zinc finger. Members of the THAP (thanatos-associated protein) family of proteins contain a well conserved DNA-binding domain known as the THAP-type zinc finger motif. Proteins containing the THAP-type zinc finger motif are commonly involved in transcriptional regulation, cell-cycle control, apoptosis and chromatin modification

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only