

Recombinant Human THAP4

Catalog No: #GP12104



Package Size: #GP12104-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human THAP4 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to C terminal 252 amino acids of human THAP4 |
| Target Name | THAP domain containing 4 |
| Other Names | PP238; CGI-36 |
| Accession No. | Swissprot:Q8WY91Gene Accession:BC000247 |
| Uniprot | Q8WY91 |
| GeneID | 51078; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

THAP domain-containing protein 4 is a protein that in humans is encoded by the THAP4 gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only