

Recombinant Human TXN2

Catalog No: #GP12106

Package Size: #GP12106-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TXN2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	thioredoxin 2
Other Names	MTRX; TRX2; MT-TRX
Accession No.	Swissprot:Q99757Gene Accession:BC013726
Uniprot	Q99757
GeneID	25828;
Storage	-20~-80°C, pH 7.6 PBS

Background

This nuclear gene encodes a mitochondrial member of the thioredoxin family, a group of small multifunctional redox-active proteins. The encoded protein may play important roles in the regulation of the mitochondrial membrane potential and in protection against oxidant-induced apoptosis.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only