

## Recombinant Human TRIM16

Catalog No: #GP12118



Package Size: #GP12118-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TRIM16
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human tripartite motif containing 16
Target Name	tripartite motif containing 16
Other Names	EBBP
Accession No.	Swissprot:O95361Gene Accession:BC053514
Uniprot	O95361
GeneID	10626;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene was identified as an estrogen and anti-estrogen regulated gene in epithelial cells stably expressing estrogen receptor. The protein encoded by this gene contains two B box domains and a coiled-coiled region that are characteristic of the B box zinc finger protein family. The proteins of this family have been reported to be involved in a variety of biological processes including cell growth, differentiation and pathogenesis. Expression of this gene was detected in most tissues. Its function, however, has not yet been determined.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only