

Recombinant Human TRIM49

Catalog No: #GP12122



Package Size: #GP12122-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TRIM49
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 200 amino acids of human tripartite motif containing 49
Target Name	tripartite motif containing 49
Other Names	RNF18; TRIM49A; TRIM49L2
Accession No.	Swissprot:P0CI25Gene Accession:BC126472
Uniprot	P0CI25
GeneID	57093;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. This gene has been found to be preferentially expressed in testis. Related pseudogenes and gene duplicates have also been identified on chromosome 11.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only