

Recombinant Human HTATIP2

Catalog No: #GP12124



Package Size: #GP12124-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human HTATIP2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 34-276 amino acids of human HIV-1 Tat interactive protein 2, 30kDa
Target Name	HIV-1 Tat interactive protein 2, 30kDa
Other Names	CC3; TIP30; SDR44U1
Accession No.	Swissprot:Q9BUP3Gene Accession:BC002439
Uniprot	Q9BUP3
GenID	10553;
Storage	-20~-80°C, pH 7.6 PBS

Background

TIP30 (30 kDa HIV-1 TAT-interacting protein), also known as HTATIP2 (HIV-1 TAT-interactive protein 2) or CC3, is a ubiquitously expressed protein that functions as an oxidoreductase. It acts as a tumor suppressor participating in metabolic suppression, inhibition of angiogenesis and induction of apoptosis.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only