

# Recombinant Human PLIN3

Catalog No: #GP12125



Package Size: #GP12125-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human PLIN3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 235 amino acids of human perilipin 3
Target Name	perilipin 3
Other Names	PP17; TIP47; M6PRBP1
Accession No.	Swissprot:O60664Gene Accession:BC005818
Uniprot	O60664
GeneID	10226;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Mannose 6-phosphate receptors (MPRs) deliver lysosomal hydrolase from the Golgi to endosomes and then return to the Golgi complex. The protein encoded by this gene interacts with the cytoplasmic domains of both cation-independent and cation-dependent MPRs, and is required for endosome-to-Golgi transport.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only