

Recombinant Human TMED3

Catalog No: #GP12132



Package Size: #GP12132-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TMED3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 194 amino acids of human transmembrane emp24 protein transport domain containing 3
Target Name	transmembrane emp24 protein transport domain containing 3
Other Names	p26; P24B; C15orf22
Accession No.	Swissprot:Q9Y3Q3Gene Accession:BC000027
Uniprot	Q9Y3Q3
GenID	23423;
Storage	-20~-80°C, pH 7.6 PBS

Background

Potential role in vesicular protein trafficking, mainly in the early secretory pathway. Contributes to the coupled localization of TMED2 and TMED10 in the cis-Golgi network.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only