

## Recombinant Human TMED4

Catalog No: #GP12133



Package Size: #GP12133-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TMED4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 198 amino acids of human transmembrane emp24 protein transport domain containing 4
Target Name	transmembrane emp24 protein transport domain containing 4
Other Names	HNLF; ERS25; GMP25iso; p24alpha3
Accession No.	Swissprot:Q7Z7H5Gene Accession:BC057851
Uniprot	Q7Z7H5
GenID	222068;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Involved in vesicular protein trafficking, mainly in the early secretory pathway. targeting. Involved in the maintenance of the Golgi apparatus. Appears to play a role in the biosynthesis of secreted cargo including processing. Involved in endoplasmic reticulum stress response. May play a role in the regulation of heat-shock response and apoptosis (By similarity).

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only