

Recombinant Human TMEM161A

Catalog No: #GP12134



Package Size: #GP12134-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TMEM161A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to N terminal 134 amino acids of human transmembrane protein 161A
Target Name	transmembrane protein 161A
Other Names	AROS-29
Accession No.	Swissprot:Q9NX61 Gene Accession:BC005210
Uniprot	Q9NX61
GeneID	54929;
Storage	-20~-80°C, pH 7.6 PBS

Background

May play a role in protection against oxidative stress. Overexpression leads to reduced levels of oxidant-induced DNA damage and apoptosis.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only