

Recombinant Human Tmprss5

Catalog No: #GP12136



Package Size: #GP12136-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Tmprss5
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to C terminal 225 amino acids of human transmembrane protease, serine 5
Target Name	transmembrane protease, serine 5
Other Names	SPINESIN
Accession No.	Swissprot:Q9H3S3Gene Accession:BC121803
Uniprot	Q9H3S3
GeneID	80975;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that belongs to the serine protease family. Serine proteases are known to be involved in many physiological and pathological processes. Alternative splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only