

Recombinant Human TMPRSS11E

Catalog No: #GP12138



Package Size: #GP12138-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TMPRSS11E
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 200-350 amino acids of human transmembrane protease, serine 11E
Target Name	transmembrane protease, serine 11E
Other Names	DESC1; TMPRSS11E2
Accession No.	Swissprot:Q9UL52Gene Accession:BC113412
Uniprot	Q9UL52
GenID	28983;
Storage	-20~-80°C, pH 7.6 PBS

Background

Serine protease which possesses both gelatinolytic and caseinolytic activities. Shows a preference for Arg in the P1 position (By similarity).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only