

## Recombinant Human TMPRSS11A

Catalog No: #GP12139



Package Size: #GP12139-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TMPRSS11A
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 203-353 amino acids of human transmembrane protease, serine 11A
Target Name	transmembrane protease, serine 11A
Other Names	ECRG1
Accession No.	Swissprot:Q6ZMR5Gene Accession:BC111796
Uniprot	Q6ZMR5
GenID	339967;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Probable serine protease which may play a role in cellular senescence. Overexpression inhibits cell growth and induce G1 cell cycle arrest.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only