

# Recombinant Human TREM1

Catalog No: #GP12156



Package Size: #GP12156-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TREM1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	triggering receptor expressed on myeloid cells 1
Other Names	CD354; TREM-1
Accession No.	Swissprot:Q9NP99Gene Accession:BC017773
Uniprot	Q9NP99
GeneID	54210;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a receptor belonging to the Ig superfamily that is expressed on myeloid cells. This protein amplifies neutrophil and monocyte-mediated inflammatory responses triggered by bacterial and fungal infections by stimulating release of pro-inflammatory chemokines and cytokines, as well as increased surface expression of cell activation markers. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only