

Recombinant Human TRIB2

Catalog No: #GP12157



Package Size: #GP12157-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human TRIB2 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Full length fusion protein |
| Target Name | tribbles pseudokinase 2 |
| Other Names | C5FW; TRB2; GS3955 |
| Accession No. | Swissprot:Q92519Gene Accession:BC002637 |
| Uniprot | Q92519 |
| GeneID | 28951; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

This gene encodes one of three members of the Tribbles family. The Tribbles members share a Trb domain, which is homologous to protein serine-threonine kinases, but lacks the active site lysine and probably lacks a catalytic function. The Tribbles proteins interact and modulate the activity of signal transduction pathways in a number of physiological and pathological processes.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only