

Recombinant Human TACSTD2

Catalog No: #GP12164



Package Size: #GP12164-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TACSTD2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 27-323 amino acids of human tumor-associated calcium signal transducer 2
Target Name	tumor-associated calcium signal transducer 2
Other Names	EGP1; GP50; M1S1; EGP-1; TROP2; GA7331; GA733-1
Accession No.	Swissprot:P09758Gene Accession:BC009409
Uniprot	P09758
GeneID	4070;
Storage	-20~-80°C, pH 7.6 PBS

Background

This intronless gene encodes a carcinoma-associated antigen. This antigen is a cell surface receptor that transduces calcium signals. Mutations of this gene have been associated with gelatinous drop-like corneal dystrophy.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only